**SureFast® Norovirus I & II 3plex (100 Reakt.)**

Art. Nr. F7140

Version 1.1

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Description

The SureFast® Norovirus I & II 3plex is a multiplex real-time RT-PCR for the direct qualitative detection and differentiation of norovirus genogroup I and genogroup II. The test contains an Internal Control RNA (ICR, consisting of MS2-bacteriophage) as an internal control of sample preparation procedure and to monitor possible PCR-inhibition. The real-time PCR assay can be used with established real-time PCR instruments, equipped for detection of three fluorescence emissions at 522 nm, 553 nm and 670 nm (FAM, VIC and Cy5) at the same time. The technical validation of instruments was performed on Agilent Mx3005P, BioRad CFX 96, Roche LightCycler® 480 II, Applied Biosystems 7500, Qiagen RotorGene Q and Cepheid SmartCycler®.

Limit of Detection

The SureFast® Norovirus I & II 3plex real-time RT-PCR has a limit of detection of ≤ 25 RNA copies. The assay limit of detection depends on RNA preparation and RNA content.

RNA-preparation

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus is recommended.

Kit components and storage

2x	Reaction Mix (700 µl)	(Code 1)
1x	PP-Mix (770 µl)	(Code 2)
1x	Enzyme Mix (80 µl)	(Code 3)
2x	Internal Control RNA (1800 µl)	(Code R)
1x	PCR Water (500 µl)	(Code N)
1x	Positive Control (100 µl)	(Code P)

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- RNA-Extraction kit (e.g. SureFast® PREP DNA/RNA Virus)
- Real-time PCR instrument, equipped with three detection channels 522 nm, 553 nm and 670 nm
- Real-time PCR consumables (plates, tubes, foils, caps)
- Pipettes with filter tips
- disposal gloves
- Vortexmixer
- Centrifuge with a rotor for the reaction tubes

Protocol**1. RNA-preparation**

For RNA-preparation the use of SureFast® PREP DNA/RNA Virus is recommended.

The test assay contains an Internal Control RNA (ICR), which can either be used as PCR inhibition control or as extraction control for the sample preparation procedure and as a PCR inhibition control.

If the ICR is used only as a PCR inhibition control, 1 µl of the ICR should be added to the Master-Mix.

If the ICR is used as an extraction control for the sample preparation procedure and as PCR inhibition control, 20 µl of the ICR should be added during extraction procedure. The ICR should always be added to the specimen-lysis buffer mixture and must not be added directly to the specimen.

2. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and blank of extraction. The test assay contains an Internal Control RNA (ICR), which can either be used as PCR inhibition control or as extraction control.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions for ICR as extraction and PCR inhibition control:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
Enzyme Mix	0.7 µl	7.7 µl
Total volume	20.1 µl	221.1 µl

Mix each master-mix well and centrifuge shortly before use.

Example for the calculation and preparation of 10 reactions for ICR only as PCR inhibition control:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	12.5 µl	137.5 µl
PP-Mix	6.9 µl	75.9 µl
Internal Control RNA	1.0 µl	11.0 µl
Enzyme Mix	0.7 µl	7.7 µl
Total volume	21.1 µl	232.1 µl

Mix each master-mix well and centrifuge shortly before use.

3. Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes or wells.
- For the negative control pipette 5 µl of PCR Water into the designated tubes and close them.
- Pipette 5 µl of sample RNA into the designated tubes and close them.
- Pipette 5 µl of the Positive Control into the designated tubes and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

4. Setup

	Rotor-Gene Q	Blockcycler
Reverse Transcription	10 min, 58°C	10 min, 58°C
Initial Denaturation (HOLD)	1 min, 95°C	1 min, 95°C
Cycles	45	45
Denaturation	15 sec, 95°C	15 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 55°C	30 sec, 55°C
Temperature Transition Rate/ Ramp Rate	Maximum	
Fluorescence Detection Setup	Detection: End of extension phase Detection system norovirus genogroup II: diverse devices FAM-channel, Quencher: BHQ Rotor-Gene Q Green LC480 II 465 nm, 510 nm Internal Control RNA: diverse devices VIC/HEX-channel, Quencher: BHQ Rotor-Gene Q Yellow LC480 II 533 nm, 580 nm Detection system norovirus genogroup I: diverse devices Cy5-channel, Quencher: BHQ Rotor-Gene Q Red LC480 II 618 nm, 660 nm	
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads/		

Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive** for norovirus genogroup II, if the sample RNA shows an amplification in the FAM-channel. A sample is stated **negative** for norovirus genogroup II, if the sample RNA shows no amplification in the FAM-channel and the internal amplification control/extraction control (VIC/HEX-channel) of the sample is **positive**.

A sample is stated **positive** for norovirus genogroup I, if the sample RNA shows an amplification in the Cy5-channel. A sample is stated **negative** for norovirus genogroup I, if the sample RNA shows no amplification in the Cy5-channel and the internal amplification control/extraction control (VIC/HEX-channel) of the sample is **positive**.

If the sample RNA and the internal amplification control are **negative** the sample contains PCR-inhibiting substances or the sample preparation was not successful. Under these circumstances an evaluation of the samples is not possible. RNA isolation and purification for the sample need to be improved.

Product Information

- Validation Report

Technical Support

For further questions please send an e-mail to info@congen.de.