RIDA[®] QUICK Aflatoxin RQS

Art. No.: R5205

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RIDA[®] QUICK Aflatoxin RQS

Brief information

RIDA[®] QUICK Aflatoxin RQS (Art. No.: R5205) is an immunochromatographic test for the quantitative determination of aflatoxin in corn. All reagents required for the assay are contained in the test kit.

Results are evaluated with the RIDA[®] QUICK SCAN reader (RQS, Art. No. ZG5005) or the RIDA[®] SMART APP software (Art. No. ZRSAM1000).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) approx. 10 min
test implementation (incubation time) 5 min

Detection limit: approx. 4 µg/kg (ppb)

Specificity: The RIDA[®] QUICK Aflatoxin RQS test determines aflatoxin in corn.

The specificity of the RIDA[®] QUICK Aflatoxin RQS test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products

RIDA [®] QUICK Aflatoxin RQS ECO	(Art. No. R5206)
RIDA [®] QUICK Aflatoxin	(Art. No. R5204)
RIDASCREEN [®] Aflatoxin Total	(Art. No. R4701)
RIDASCREEN [®] Aflatoxin B1 30/15	(Art. No. R1211)
RIDASCREEN [®] FAST Aflatoxin	(Art. No. R5202)
RIDASCREEN [®] FAST Aflatoxin SC	(Art. No. R9002)
RIDASCREEN [®] Aflatoxin M1	(Art. No. R1121)
RIDASCREEN [®] FAST Aflatoxin M1	(Art. No. R5812)

1. Intended use

RIDA®QUICK Aflatoxin RQS is a quantitative immunochromatographic test in strip format for the determination of aflatoxin in corn.

2. General

Aflatoxins are secondary metabolites of the fungi species *Aspergillus flavus*, *parasiticus* and *nomius*. These fungi occur in humid tropical areas and the contamination of agricultural products takes place in the cultivable countries. Aflatoxins belong to the strongest natural occurring cancerogenic substances.

3. Test principle

The basis of the immunochromatographic assay in test strip formate is an antigen-antibody reaction. A specific antibody against aflatoxin recognizes the aflatoxin molecules in the samples.

The control band (control line) should be present in all cases in order to prove that the test strip is valid. The control band intensity decreases with increasing aflatoxin concentration in the sample. When a sample contains very high levels of aflatoxin, there is only a very weak control band visible.

The intensity of the test band (test line) depends on the aflatoxin concentration of the sample.

Results are evaluated with the RIDA®QUICK SCAN reader or the RIDA®SMART APP software.

4. Reagents provided

Each kit contains sufficient materials for 20 determinations. Each test kit contains:

Component	Cap color	Format	Content
20 x Test strip	-	Ready to use	One determination per strip (separately packed)
1 x Mobile solvent	White	Ready to use	4.5 ml
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for use with RIDA®SMART APP

5. Materials required but not provided

5.1. Equipment:

- laboratory mincer / grinder
- balance
- optional: shaker
- optional: filter paper, Whatman No. 1 (equivalent or centrifuge)
- graduated cylinder
- 50 and 100 µl pipette
- (e.g. R-Biopharm FP 50 / FP 100, Art. No.: Z0006 / Z0007)
- RIDA[®]QUICK SCAN (Art. No.: ZG5005) or activated RIDA[®]SMART APP software (Art. No. ZRSAM1000) installed on compatible smartphone

5.2. Reagents:

- 70 % methanol (e.g. 70 ml methanol (100%) / 30 ml distilled water), exact adjustment necessary
- 50 % ethanol (e.g. 50 ml ethanol (100 %) / 50 ml distilled water), exact adjustment necessary

6. Warnings and precautions for the users

This test should only be carried out by trained employees. The instruction for use must be strictly followed.

The test strips are sensitive to humidity. Humid test strips may influence the test results negatively. For this reason keep the strips away from humidity. This has to be noted especially for already opened test strip packings.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at (2 - 8 °C / 35 - 46 °F). Do not freeze components of the test kit.

No quality guarantee is accepted after expiry of the kit on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Sample preparation

Bring all components to room temperature (20 - 25 °C / 68 - 77 °F) before use!

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

8.1. USDA / GIPSA extraction method with ethanol

(Grain Inspection, Packers and Stockyards Administration (GIPSA))

- weigh 50 g of ground sample into a suitable vial
- transfer 50 g of this sample into a clean blender jar
- add 100 ml 50 % ethanol
- cover the blender jar and blend for 2 min (maximum speed)
- allow the sample to settle for 1 min
- transfer 1.0 ml of extraction mixture to a centrifuge tube
- centrifuge for 1 min
- pipette 200 µl of mobile solvent and 100 µl of the centrifuged extract into a clean tube (use a new pipette tip for each sample) and close the tube
- invert tube 5 times to mix thoroughly
- use 100 µl in the test

8.2. Sample extraction with methanol or ethanol

- weigh 10 g of ground sample into a suitable vial and add 20 ml 70 % methanol or ml 50 % Ethanol
- close the tube and shake the sample vigorously for 5 min (manually or with shaker / vortex)
- let the solution come to sedimentation for 3 min, filtrate or centrifuge (it is absolutely necessary to avoid methanol or ethanol evaporation)
- mix 100 µl of the temperate mobile solvent with 50 µl of the clear supernatant and use 100 µl in the test

Remark:

If necessary the sample weight can be changed. Then the volume of methanol / ethanol has to be adapted, e.g. 5 g in 10 ml methanol (70 %) or 10 ml ethanol (50 %).

9. Test procedure

- apply 100 µl of the sample solution (see 8. Sample preparation) on the application area of the test strip
- evaluate the result after 5 min incubation with RIDA[®]SMART APP or RIDA[®]QUICK SCAN

10. Evaluation


The right band in the test strip reaction field is a control band (control line) and must appear after each test procedure (see Fig. 1). If the band is missing, the test result is not valid because of improper test procedure or deterioration of the reagents. Repeat the test with a new strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is contaminated > 150 µg/kg visibility of the control band is very weak and the test may give an invalid result.

The left band in the reaction field is the test band (test line). Its intensity depends on the aflatoxin concentration of the sample.

To obtain comparable results, the test strips must always be evaluated after 5 min.

Fig.1: Test strip RIDA[®]QUICK Aflatoxin RQS

10.1. Evaluation with RIDA®QUICK SCAN (Art. No. ZG5005)

First, please read the instructions for use of the RIDA®QUICK SCAN reader attentively. The test strip evaluation has to be performed according to these guidelines.

Please find a short description of how to use of the RIDA®QUICK SCAN reader in the test kit.

Please note:

The test kit contains several quality assurance certificates. Please choose the corresponding certificate for the used extraction method. Please use the parameters on the appropriate certificate (barcode with check sum) for the analysis with the RIDA®QUICK SCAN.

10.2. Evaluation with RIDA®SMART APP (Art. No. ZRSAM1000)

First, please read the instructions for use of the RIDA®SMART APP software attentively. The test strip evaluation has to be performed according to these guidelines.

Please find a short description of how to use the RIDA®SMART APP on our website www.r-biopharm.de.

For the test strip evaluation the RIDA®SMART APP cover is needed which is included in the test kit. The RIDA®SMART APP cover contains lot-specific information and must be used only with the appropriate lot of the test kit. For further information please check the certificate of analysis enclosed in the kit.

Before starting the test implementation and test evaluation please check if the correct RIDA®SMART APP cover is available.

RIDA®SMART APP software applications for the evaluation of RIDA®QUICK Aflatoxin RQS test strips:

Extraction method	RIDA®SMART APP application
70 % methanol	Methanol
50 % ethanol	Ethanol

11. Sensitivity

The RIDA[®]QUICK Aflatoxin RQS test is capable of detecting aflatoxin contaminations of $\geq 4 - 100 \mu\text{g}/\text{kg}$.

For further product information and applications please contact your local distributor or R-Biopharm AG (sales@r-biopharm.de).

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG