# **RIDA<sup>®</sup> QUICK Fumonisin RQS**

**Art. No.: R5606**

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# RIDA<sup>®</sup> QUICK Fumonisin RQS

## Brief information

RIDA<sup>®</sup> QUICK Fumonisin RQS (Art. No.: R5606) is an immunochromatographic test for the determination of fumonisin in corn.

All reagents required for the assay are contained in the test kit.

Results are evaluated with with the RIDA<sup>®</sup> QUICK SCAN (RQS, Art. No. ZG5005) or the RIDA<sup>®</sup> SMART APP software (Art. No. ZRSAM1000).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) .....approx. 10 min  
test implementation (incubation time) .....5 min

Limit of detection 0.3 mg/kg (ppm) - see 8.1.  
3.0 mg/kg (ppm) - see 8.2.

Specificity: The RIDA<sup>®</sup> QUICK Fumonisin test determines fumonisin in corn samples. The test has been validated for this matrix.

The specificity of the RIDA<sup>®</sup> QUICK Fumonisin RQS test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

## Related products

RIDASCREEN<sup>®</sup> Fumonisin (Art. No. R3401)

RIDASCREEN<sup>®</sup> FAST Fumonisin (Art. No. R5602)

## 1. Intended use

RIDA®QUICK Fumonisin RQS is a quantitative immunochromatographic test in strip format for the determination of fumonisin in corn.

## 2. General

Fumonisin is carcinogenic, highly toxic metabolite of *Fusarium moniliforme*, which grows host-specific on corn. The toxic concentrations of fumonisin differ significantly depending on the animal species.

## 3. Test principle

The basis of the immunochromatographic assay in test strip format is an antigen-antibody reaction. A specific antibody against fumonisin binds the fumonisin molecules in the samples.

Results are evaluated with RIDA®QUICK SCAN (RQS) reader or the RIDA®SMART APP software.

The control band (control line) should be present in all cases in order to prove that the test strip is valid.

## 4. Reagents provided

Each kit contains sufficient materials for 20 determinations.

Each test kit contains:

Component	Cap color	Format	Content
20 x Test strip	-	Ready to use	One determination each (separately packed)
1 x Extraction buffer	Transparent	Concentrate 5 x	100 ml
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for the use with RIDA®SMART APP

## **5. Materials required but not provided**

### 5.1. Equipment:

- laboratory mincer / grinder
- balance
- optional: shaker
- optional: filter paper (Whatman No. 1), equivalent or centrifuge
- graduated cylinder
- 100 µl, 300 µl and 700 µl micropipettes  
(e.g. R-Biopharm FP 100, Art. No.: Z0007)
- RIDA<sup>®</sup>QUICK SCAN (Art. No.: ZG5005) or activated RIDA<sup>®</sup>SMART APP software (Art. No. ZRSAM1000) installed on compatible smartphone

## **6. Warnings and precautions for the users**

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The test strips are sensitive to humidity. Humid test strips may influence the test results negatively; therefore keep the strips away from humidity. This is especially important for already opened test strip packings.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## **7. Storage instructions**

Store the kit at 2 - 8 °C / 36 - 47 °F. DO NOT FREEZE the test strips.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

## 8. Sample preparation

Dilute the extraction buffer 1:5 (1+4) with distilled water. Store this ready to use buffer at 2 – 8 °C / 35 – 46 °F. The expiry is not affected by dilution.

**Bring test strips and ready to use extraction buffer to room temperature (20 - 25 °C / 68 - 77 °F) before use!**

The samples should be stored in a cool place, protected from light.

A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

### 8.1. Detection limit 0.3 mg/kg

- weigh 5 g of ground sample into a screw cap tube and add 20 ml of ready to use extraction buffer
- close the tube and shake the sample vigorously for approx. 5 min (manually or with shaker)
- let the solution settle for approx. 3 - 5 min, filtrate or centrifuge
- mix 100 µl of the clear supernatant with 700 µl ready to use extraction buffer (= dilution I)
- apply 100 µl of dilution I onto the application area of the test strip

### 8.2. Detection limit 3.0 mg/kg

- weigh 5 g of ground sample into a screw cap tube and add 20 ml of ready to use extraction buffer
- close the tube and shake the sample vigorously for approx. 5 min (manually or with shaker)
- let the solution settle for approx. 3 - 5 min, filtrate or centrifuge
- mix 100 µl of the clear supernatant with 700 µl ready to use extraction buffer (= dilution I)
- mix 100 µl of dilution I with 300 µl ready to use extraction buffer (= dilution II)
- apply 100 µl of dilution II onto the application area of the test strip

## 9. Test procedure

- apply 100 µl of dilution I or dilution II (see 8.1. or 8.2.) onto the application area of the test strip (see Fig. 1)
- read the result after 5 min (+ max. 15 sec) incubation with RIDA<sup>®</sup>QUICK SCAN or RIDA<sup>®</sup>SMART APP


## 10. Evaluation

The right band in the test strip reaction field is a control band (control line) and must appear after each test procedure. If the band is missing, the test result is not valid because of improper test procedure or deterioration of the reagents. Repeat the test with a new strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is contaminated > 25 mg/kg visibility of the control band is very weak.

To obtain comparable results, the test strips must always be evaluated after 5 min.

Fig.1: Test strip RIDA<sup>®</sup>QUICK Fumonisin RQS

### 10.1. Evaluation with RIDA<sup>®</sup>QUICK SCAN (Art. No. ZG5005)

First, please read the instructions for use of the RIDA<sup>®</sup>QUICK SCAN reader attentively. The test strip evaluation has to be performed according to these guidelines.

RIDA<sup>®</sup>QUICK SCAN software methods for the evaluation of RIDA<sup>®</sup>QUICK Fumonisin RQS test strips:

Detection range	RIDA <sup>®</sup> QUICK SCAN method
0.3 – 3.0 mg/kg	FUMO 0.3-3.0
3.0 – 10.0 mg/kg	FUMO 3.0-10.0

## 10.2. Evaluation with RIDA<sup>®</sup>SMART APP (Art. No. ZRSAM1000)

First, please read the instructions for use of the RIDA<sup>®</sup>SMART APP software attentively. The test strip evaluation has to be performed according to these guidelines.

Please find a short description of how to use the RIDA<sup>®</sup>SMART APP on our website [www.r-biopharm.de](http://www.r-biopharm.de).

For the test strip evaluation the RIDA<sup>®</sup>SMART APP cover is needed which is included in the test kit. The RIDA<sup>®</sup>SMART APP cover contains lot-specific information and must be used only with the appropriate lot of the R5606 RIDA<sup>®</sup>QUICK Fumonisin RQS test kit. For further information please check the certificate of analysis enclosed in the kit.

Before starting the test implementation and test evaluation please check if the correct RIDA<sup>®</sup>SMART APP cover is available.

RIDA<sup>®</sup>SMART APP software applications for the evaluation of RIDA<sup>®</sup>QUICK Fumonisin RQS test strips:

Detection range	RIDA <sup>®</sup> QUICK SCAN application
0.3 – 3.0 mg/kg	0,3-3,0ppm
3.0 – 10.0 mg/kg	3,0-10,0ppm

## 11. Sensitivity

The RIDA<sup>®</sup>QUICK Fumonisin RQS test is able to detect fumonisin contaminations of  $\geq 0.3$  mg/kg to 3 mg/kg and of 3 mg/kg to 10 mg/kg.

**For further product information please contact your local distributor or R-Biopharm AG ([sales@r-biopharm.de](mailto:sales@r-biopharm.de)).**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG