# **RIDA<sup>®</sup> QUICK Gliadin (single packaged)**

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Immuno chromatographic test for the detection of Gliadin and corresponding prolamins on surfaces and in food

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**Approved as AOAC Official Method of Analysis (OMA)**

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Please contact for questions and further information:

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RIDA<sup>®</sup> and RIDASCREEN<sup>®</sup>  
are registered trademarks of R-Biopharm AG  
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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# RIDA<sup>®</sup>QUICK Gliadin (single packaged)

## Brief information

RIDA<sup>®</sup>QUICK Gliadin (single packaged) (Art. No. R7004) is an immunochromatographic test for the qualitative detection of gliadin / gluten contamination.

- on surfaces (swab test for the hygiene control in production and in laboratories)
- in gluten-free raw material after an ethanol extraction
- in gluten-free processed food after an extraction with the Cocktail (patented).

**The R5 dip stick RIDA<sup>®</sup>QUICK Gliadin has been approved as AOAC-OMA.**

All reagents required for the swab test are contained in the test kit.

The test kit contains 25 test strips (single packaged) for 1 determination each. Results are evaluated visually.

Time requirement:      sampling for swab test..... approx. 1 min  
                                 sample preparation  
                                 for 10 raw materials ..... approx. 15 min  
                                 for processed food..... approx. 120 min  
                                 test implementation (incubation time) ..... 5 min

Detection limit:      - on **surfaces** approx. 1 - 2 µg gliadin / 100 cm<sup>2</sup>  
                                 (approx. 2 - 4 µg gluten / 100 cm<sup>2</sup>)  
                                 - in **raw material** approx. 2.2 mg/kg gliadin  
                                 (approx. 4.4 mg/kg gluten) [Prediction Interval 3.5 - 5.6]  
                                 - in **processed food** approx. 3.1 mg/kg gliadin  
                                 (approx. 6.3 mg/kg gluten) [Prediction Interval 3.9 - 10.2]

Specificity:              The **monoclonal antibody R5** reacts with the gliadin-fraction from wheat and corresponding prolamins from rye and barley.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

## Related products:

RIDASCREEN® Gliadin (Art. No. R7001)  
RIDASCREEN® FAST Gliadin (Art. No. R7002)  
RIDASCREEN® Gliadin competitive (Art. No. R7021)  
RIDA® QUICK Gliadin (Art. No. R7003)  
RIDA® QUICK Gliadin (ready to swab) (Art. No. R7005)  
Cocktail (patented) (Art. No. R7006 / R7016)  
RIDA® Extraction Solution (colorless) (Art. No. R7098)  
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)  
SureFood® ALLERGEN real time PCR Gluten (Art. No. S3106)  
SureFood® ALLERGEN QUANT real time PCR Gluten (Art. No. S3206)

## 1. Intended use

RIDA® QUICK Gliadin can be used as a swab test for the gluten determination on surfaces in the hygiene control and for the qualitative detection of gliadin / gluten in raw material and processed food. The test has been developed for the detection of low amounts of gluten (contamination). **No** high-dose-hook-effect is observed at high concentrations. However, the red target band may smear at high gluten concentrations.

## 2. General

The use of wheat flour and gluten in foodstuff is extremely common because of their useful effects on e.g. texture, moisture retention and flavour. Gluten is a mixture of prolamin and glutelin proteins present in wheat, rye and barley.

Coeliac disease is a permanent intolerance to gluten that results in damage to the small intestine and is reversible when gluten is avoided by diet.

The Codex Alimentarius Commission has stipulated in the „Codex Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten” (CODEX STAN 118-1979) the limit value for gluten-free food at 20 mg/kg gluten. This threshold has also been adopted by many national legislations. The prolamin content (e.g. gliadin) of gluten is generally taken as 50 % (CODEX STAN 118-1979).

The official type I method for gluten determination according to the Codex Alimentarius is an ELISA which uses the R5 antibody (Mendez). This requirement is fulfilled by RIDASCREEN® Gliadin test (Art. Nr. R7001). **The RIDA® QUICK Gliadin test strips also use the R5 antibody and show a good correlation with the official method, the R5-ELISA RIDASCREEN® Gliadin. R-Biopharm AG is the only company, that is allowed to use the R5 antibody for the test strips.**

### 3. Test principle

The basis of the immunochromatographic test is the monoclonal R5-antibody which is specific for the detection of gliadin from wheat and prolamins from rye and barley. If gliadin is present, a sandwich is formed consisting of immobilized R5 antibody at the target band, gliadin and red latex-labeled R5 antibody. Results are read visually. Generally, the higher the analyte level in the sample, the stronger the red color of the test band will be.

### 4. Reagents provided

Each kit contains sufficient materials for 25 measurements.

Component	Cap color	Format	Volume
Test strip	-	Ready to use, single packaged	25 pieces
Test tube			30 pieces
Buffer	transparent	Ready to use	60 ml
Evaluation card			1 piece

### 5. Materials required but not provided

#### 5.1. Equipment:

##### **For analysis of raw material and processed food**

- scales
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- shaker
- centrifugal vials + centrifuge or paper filter
- graduated pipettes

#### 5.2. Reagents:

##### **For analysis of raw material**

- distilled or deionized water
- ethanol solution (60 %), for the extraction of the samples (add 150 ml ethanol p.a. to 100 ml distilled water and shake well)
- for soy containing food: add 1 g skim milk powder (food quality) to 1 g sample

## **For analysis of processed food**

- distilled or deionized water
- ethanol solution (80 %), for the extraction of the samples (add 120 ml ethanol p.a. to 30 ml distilled water and shake well)
- Cocktail (patented) (Art. No. R7006 / R7016)
- for polyphenol containing food: add 0.25 g skim milk powder (food quality) to 0.25 g sample

## **6. Warnings and precautions for the users**

Airborne cereal dust and dirty laboratory equipment lead to gliadin contamination of the assay. In order to avoid cross-contamination by cereal dust, please note the following points:

- wear gloves before starting and during the assay
- clean surfaces, glass vials, mincers and other equipment with 40 % ethanol or 2-propanol
- for the analysis of raw material, the extraction and the test procedure should be carried out in separate rooms

The dip sticks are very sensitive to humidity that could turn the test useless. For this reason keep the strips away from humidity!

This kit may contain further hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## **7. Storage instructions**

Store the unopened kit at 2 - 8 °C (36 - 46 °F). Do not freeze the kit.

No quality guarantee is accepted after the expiry date on the kit label.

## **8. Test procedure**

### **8.1. Swab test: sampling and test implementation**

1. Take as many test tubes as samples to be analyzed.
2. Place 500 µl of buffer in the test tube (e.g. using the disposable pipette provided).

3. Swab the lower end (reaction zone) of a dry dip stick thoroughly over a sampling area of 10 x 10 cm (wear gloves).
4. Place the dip stick vertically into the test tube with the arrow pointing down. Do not immerse the dip stick beyond the maximum line (look at figure under 8.4).
5. Take out the strip after exactly 5 min (+/- 10 s) and read the result using the evaluation card.

## 8.2. Sample preparation raw material (non-processed food)

### 8.2.1. Fluid and soft raw material

- **fluid raw material:** mix 1 ml of the sample with 9 ml 60 % ethanol solution
- for soy milk add additionally 1 g of skim milk powder
- **soft raw material:** weigh 1 g of a representative sample and add 10 ml 60 % ethanol solution
- for soy milk products add additionally 1 g of skim milk powder
- shake well for at least 30 sec. (vortex)
- centrifuge: 10 min / at least 2500 g / room temperature (20 - 25 °C / 68 - 77 °F)
- alternative: let the sample settle down and / or filtrate

### 8.2.2. Solid and hard raw material

- weigh 5 g sample and grind it to powder
- use 1 g of this powder and add 10 ml 60 % ethanol solution (for soy containing samples add 1 g of skim milk powder)
- shake well for at least 30 sec. (vortex)
- centrifuge: 10 min / at least 2500 g / room temperature (20 - 25 °C / 68 - 77 °F)
- alternative: let the sample settle down and / or filtrate

### 8.3. Sample preparation processed food with Cocktail (patented)

Homogenize well a sufficient amount (at least 50 g or 50 ml) of sample (grind it thoroughly to powder and mix well or mix well the solution respectively).

- **liquid food samples:** use 0.25 ml of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **other food samples (e.g. soy and quinoa containing samples):** weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **tannin and polyphenol containing food samples (e.g. chocolate, coffee, cocoa, chestnut flour, buckwheat, millet and spices):** weigh 0.25 g of the homogenized sample, add 0.25 g of skimmed milk powder and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **meat and sausages:** in these matrices gliadin may be not distributed evenly, therefore, weigh 50 g sample and homogenize: weigh 0.25 g of the homogenized sample and add 2.5 ml of the Cocktail (patented), close the vial and mix well
- **oat samples:** gliadin may be not distributed evenly; furthermore the samples are difficult to homogenize. Therefore, homogenize 200g, then carry out the extraction with at least the fourfold amount of reagents: weigh 1 g of the homogenized sample and add 10 ml of the Cocktail (patented), close the vial and mix well

#### **Please further extract all samples as described in the following:**

- incubate for 40 min at 50 °C (122 °F) in the water bath
- let the sample cool down and then mix it with 7.5 ml 80 % ethanol (see 5.2.)  
(for oat samples: 30 ml 80% ethanol)
- close the vial and shake for 1 h up side down or by a rotator at room temperature (20 - 25 °C / 68 - 77 °F)
- centrifuge: 10 min, at least 2500 g, at room temperature (20 - 25 °C / 68 - 77 °F) and / or filter the extract (alternatively 2 ml of the extract can be centrifuged with high speed for 10 min in reaction caps by using a microcentrifuge)
- put the supernatant in a screw top vial

**All supernatants / filtrates obtained after centrifugation or filtration can be stored in a tightly closed vial in the dark at room temperature (20 - 25 °C / 68 - 77 °C) up to four weeks.**

#### 8.4. Test implementation for raw material and processed food

1. Take as many test tubes as samples to be analyzed.
2. Place 500 µl of buffer in the test tube (e.g. using the disposable pipette provided).
3. Pipette 50 µl of the sample supernatant / filtrate or place 3 drops with the provided disposable pipette, vertically dropped in the test tube and shake slightly.
4. Place the dip stick vertically into the test tube with the arrow pointing down. Do not immerse the strip beyond the maximum line.
5. Take out the dip stick after exactly 5 min (+/- 10 s) and read the result using the evaluation card.

## 9. Results

### **Positive result: two colored bands**

The sample is positive if two colored bands (the blue control band and the red test band) are visible within the result window. In the case of swabbing test bands with non-uniform intensity may occur due to an inhomogeneous gluten distribution on the surface or different swabbing procedures.

Swab test: > approx. 1 - 2 µg gliadin / 100 cm<sup>2</sup>  
(approx. 2 - 4 µg gluten / 100 cm<sup>2</sup>)

Raw material: > approx. 2.2 mg/kg gliadin (approx. 4.4 mg/kg gluten)

Processed food: > approx. 3.1 mg/kg gliadin (approx. 6.3 mg/kg gluten)

### **Negative result: only the blue control band**

The sample is negative if no red test band is visible within the result window.

Swab test: < approx. 1 - 2 µg gliadin / 100 cm<sup>2</sup>  
(approx. 2 - 4 µg gluten / 100 cm<sup>2</sup>)

Raw material: < approx. 2.2 mg/kg gliadin (approx. 4.4 mg/kg gluten)

Processed food: < approx. 3.1 mg/kg gliadin (approx. 6.3 mg/kg gluten)

### **Invalid result: no colored band**

If no band is visible within the result window after performing the test, the test is considered invalid.

## In general

- The test strip has been developed for the detection of gluten contamination.
- The limit of detection is dependent on sample type and extraction efficiency or the properties of the swabbed surface and the kind of contamination respectively.
- The sample extraction with ethanol should only be used for raw material that were surely not heated and not processed.
- A negative result does not necessarily indicate the absence of gluten as gluten may be not homogenously distributed or the level of gluten in the product is below the limit of detection.

## Recommendation

- For documentation, the upper part of the dip stick marked with “Gluten” together with the test bands must be cut off.
- The use of assay test controls (R7012, for cocktail extraction) or of spiked samples is recommended for quality control.
- It is recommended comparing the extraction efficiency of ethanol with the Cocktail (patented) (R7006).
- The RIDASCREEN® Gliadin (Art. No. R7001) should be used for quantification. This test kit is also AOAC-RI and AOAC-OMA (Official Method of Analysis) validated.

**Further information, validation data and applications are available on request from your local distributor or R-Biopharm AG.**

## Further applications:

- Sample preparation for processed food with the RIDA® Extraction Solution (colorless) (Art. No. R7098) - **only after validation**
- Sample preparation for polyphenol containing raw materials (e.g. chocolate, coffee, cacao, buckwheat) with fish gelatine.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

# RIDA® QUICK Gliadin (single packaged)

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RIDASCREEN® Gliadin (Art. No. R7001)  
RIDASCREEN®FAST Gliadin (Art. No. R7002)  
RIDASCREEN® Gliadin competitive (Art. No. R7021)  
RIDA®QUICK Gliadin (Art. No. R7003)  
RIDA®QUICK Gliadin (ready to swab) (Art. Nr. R7005)  
Cocktail (patented) (Art. No. R7006 / R7016)  
RIDA® Extraction Solution (colorless) (Art. No. R7098)  
Set of 3 processed Gliadin Assay Controls (Art. No. R7012)  
SureFood® ALLERGEN real time PCR Gluten (Art. No. S3106)  
SureFood® ALLERGEN QUANT real time PCR Gluten (Art. No. S3206)

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