**SureFood® ANIMAL QUANT Beef (2 x 50 Reakt.)**

Art. Nr. S1010

Version 2.3

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	Blockcycler/LightCycler® 480	Rotorcyclcr
Initial Denaturation (HOLD) Cycles Denaturation Annealing Extension (CYCLE)	5 min, 95°C 45 10 sec, 95°C 15 sec, 62°C 30 sec, 65°C	1 min, 95°C 45 5 sec, 95°C 10 sec, 62°C 15 sec, 65°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: TAMRA Passive Reference: none	LightCycler® Channel: 530/610 oder F1/F2 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detaillierte Informationen zur Einstellung bestimmter real-time PCR Geräte stehen auf der CONGEN-Homepage zur Verfügung: http://www.congen.de/unternehmen/download		

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Description

The test detects the relative beef DNA amount to the total animal DNA content in meat samples. Therefore the kit contains two PCR systems, one for detection of a beef specific gene (**Bos** detection gene) and one for the detection of an animal gene (**Ref** reference gene). The real-time PCR assay can be used with established real-time PCR instruments (Roche LightCycler®, Rotor-Gene Q, ABI PRISM, Eppendorf realplex, BioRad CFX96, Agilent MxSeries etc.).

Limit of Detection

The beef specific PCR has a limit of detection of ≤ 5 DNA-copies. The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content. The limit of quantitation depends on the concentration of the sample DNA used in the analysis. For example, if 50,000 target-sequence copies of the reference gene are present, the relative quantitation limit for beef DNA is 0.04 %.

Cross reactivity

Cross reactivity was observed with DNA extracts from *Bison sp.* (100 %) and Water buffalo (*Bubalis arnee*) with 0.6%.

DNA-preparation

For DNA-preparation of meat samples the use of SureFood® PREP Basic is recommended.

Kit components and storage

Kit Code	Reagent	Amount	Lid Color
1	Ref Reaction Mix	2 x 450 µl	Yellow
2	Bos Reaction Mix	2 x 450 µl	Yellow
3	Taq Polymerase	1 x 15 µl	Red
4	Dilution Buffer	1 x 1300 µl	White
5	Standard DNA	1 x 25 µl	Dark Blue
6	Positive Control (100% beef)	1 x 25 µl	Light Blue

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- real-time PCR instrument
- real-time PCR consumables (plates, tubes, foils, capillaries, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples, controls and standards).

Reactions needed for the reference-gene detection:

For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)

For each sample: at least 1 reaction with each sample DNA

Reactions needed for the beef detection:

For each run: 5 reactions for the standard curve
3 reactions for controls (1x no-template and 2x Positive Control)

For each sample: at least 1 reaction with each sample DNA

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Ref Reaction Mix or Bos Reaction Mix	18.0 µl	198.0 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	18.1 µl	199.1 µl

Mix each master-mix well and centrifuge shortly before use.

2. Setup

	Blockcycler/LightCycler® 480	Rotorcyclcr
Initial Denaturation (HOLD) Cycles Denaturation Annealing Extension (CYCLE)	5 min, 95°C 45 10 sec, 95°C 15 sec, 62°C 30 sec, 65°C	1 min, 95°C 45 5 sec, 95°C 10 sec, 62°C 15 sec, 65°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: TAMRA Passive Reference: none	LightCycler® Channel: 530/610 or F1/F2 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: FAM (Green)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

3. Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different DNA concentrations for the standard curves of the **reference** gene and the **detection** gene. Prepare 5 dilutions of the supplied Standard DNA with the supplied Dilution Buffer. Prepare 5 reaction tubes (labeled S1 to S5) and add 45 µl Dilution Buffer (**Code 4**) each. The following procedure is recommended:

standard	dilution	copy number per µl	final copy number per reaction*
S1	45 µl Dilution Buffer + 5 µl Standard DNA	100,000 copies	200,000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	10,000 copies	20,000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	1000 copies	2000 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	100 copies	200 copies
S5	45 µl Dilution Buffer + 5 µl DNA of S4	10 copies	20 copies

***Note:** 2 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

If the diluted DNA standards (S1 to S5) are not immediately used, store them at -20°C. The dilutions are stable up to two months at -20°C. Before use allow the reagents to thaw, mix them on a vortex and centrifuge carefully before opening and use.

4. Preparation of the PCR-mix

- Pipette 18 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 2 µl of sample DNA into the designated tubes/wells or capillaries.
- Pipette 2 µl of the Positive Control and the standard dilutions into the designated tubes/wells or capillaries.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates or capillaries into the PCR instrument and start the run according to the setup.

Interpretation of results

The calculation for both reaction systems (**Ref**, **Bos**) has to be made separately. Mark the standards, the controls and the samples for the specific system (**Bos**) and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. Repeat the same procedure for the reference gene system (**Ref**). The value for the slope of the standard curve has to be between -3.1 and -3.6 and the correlation coefficient $R^2 > 0.98$. In case of different values for the standard curve, it should not be used for calculation.

By using the calculated copy numbers for **Ref** and **Bos** the relative beef content of the sample DNA and the Positive Control can be determined in the following way (example):

Sample beef (Bos)	1350 copies	Positive Control beef (Bos)	900 copies
Sample animal (Ref)	45,000 copies	Positive Control animal (Ref)	1100 copies

Divide the copy number of the specific system by the copy number of the reference gene system and multiply by 100 to obtain the percentage.

beef DNA content = **Bos** copy number * 100 / **Ref** copy number

sample DNA **beef** content = $1350 * 100 / 45,000$ sample-DNA **beef** content = 3 %

For the given example the numbers lead to a **beef** DNA content of 3.0 %, for the sample and 81.8 % for the Positive Control with the same calculation.

For a final calculation the use of a correction factor K for the correction of run-to-run fluctuations is necessary. The correction factor is the relation of the true percentage value of the Positive Control (100 % **beef** content) and the measured percentage of the Positive Control. The factor is calculated in the following way:

$K = \text{true percentage of Positive Control} / \text{measured percentage of Positive Control}$

$K (\text{example}) = 100 \% / 81.8 \% = 1.2$

The calculated value for the sample is multiplied with K to obtain a corrected value.

percentage sample = measured percentage of sample * K sample (example) = $3.0 \% * 1.2 = 3.6 \%$

For that example the **beef** DNA content is 3.6 %.

Note

Innards, particularly kidney or liver, show a higher DNA content per weight unit when compared to muscle tissue. This could lead to over-estimation of the species content of the animal from which the innards come from. For meat products with innards content >5 % care has to be taken when relative species contents are reported.

Product Information

- Validation Report
- Microsoft Excel template of calculation
(Download: www.congen.de/en/company/downloads)

Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

Distribution and ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com