User Manual

SureFood® PREP Basic

Art. No. S1052

100 extractions

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Description

This kit is intended to be used for the extraction of animal and plant DNA (deoxyribonucleic acid) from raw materials as well as slightly processed food and feed. It is also recommended for the extraction of animal DNA from highly processed food and feed.

Principle

1. Preparation of the basic material
2. Lysis of the basic material
3. Pre-filtration and setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

Kit components (per box) and storage

1x Proteinase K	(Code K)	1x Elution Buffer (10 ml)	(Code E)
1x PCR grade water (1.1 ml)	(Code H)	1x Spin Filter, green (50 x)	(Code F)
1x Lysis Buffer (20 ml)	(Code L)	1x Spin Filter, clear (50 x)	(Code S)
1x Binding Buffer (10 ml)	(Code B)	2x Receiver Tubes 2.0 ml, clear (50 x)	(Code R)
1x Pre-Wash Buffer (60 ml)*	(Code P)	1x Receiver Tubes 1.5 ml, clear (50 x)	(Code T)
1x Wash Buffer (60 ml)*	(Code W)		

* After adding ethanol (purity \geq 96 %; not supplied with the kit)

The components are equally divided into two boxes.

The Proteinase K should be stored at -20°C. All other reagents of the kit should be stored dry and at room temperature (14-25°C).

Before every extraction make sure that all components have room temperature. If there are any precipitates within the provided solutions, solve them by warming carefully (up to 30°C).

Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 65°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol (purity \geq 96 %) for preparation of Pre-Wash Buffer and Wash Buffer
- disposal bags or waste bin

A flow chart for visualization of the extraction protocol is available at the CONGEN homepage: <http://www.congen.de/en/company/downloads>.

Preparations

General

- Add 1 ml PCR grade water (**Code H**) to the Proteinase K (**Code K**) and mix thoroughly.
- Add 30 ml ethanol to the Pre-Wash Buffer (**Code P**) and mix thoroughly.
- Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

Before each preparation

- Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer (**Code E**) under calculation of a reserve volume into a reaction tube (not provided with the kit) and equilibrate to 65°C (The Elution Buffer is necessary for step 7).

Extraction protocol

1. Preparation of the basic material

Transfer 50 mg homogenized sample material into a 2.0 ml reaction tube (not provided with the kit). For samples with high water content it is necessary to take up to 100 mg of the basic material.

2. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**) and 20 µl of Proteinase K (**Code K**) to the reaction tube and mix it briefly. Incubate on a heating block under continuously shaking for 30 minutes at 65°C.

Note: For strongly swelling samples it can be necessary to add additional volume of Lysis Buffer during the lysis.

3. Pre-filtration and setting of optimal binding conditions

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Place a green Spin Filter (**Code F – green, for pre-filtration**) into a clear 2.0 ml Receiver Tube (**Code R**). Transfer the liquid supernatant directly onto the Spin Filter.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the Spin Filter.

Note: For non-processed animal material the usage of the green Spin Filter is superfluous. Transfer the liquid supernatant directly into a clear 2.0 ml Receiver Tube (**Code R**).

4. Binding of the nucleic acids on a Spin Filter

Add 200 µl Binding Buffer (**Code B**) to the filtrate and mix.

Place a clear Spin Filter (**Code S – clear, for binding**) into a new, clear 2.0 ml Receiver Tube (**Code R**). Transfer the complete solution onto the Spin Filter. Incubate at room temperature for 1 min.

Centrifuge the Spin Filter with the Receiver Tube for 1 min at 12,000 rpm. After centrifugation discard the filtrate and place the Spin Filter back into the Receiver Tube.

5. Purification of the bound nucleic acids

Add 550 µl Pre-Wash Buffer (**Code P**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Note: For non-processed animal material the usage of the Pre-Wash Buffer is superfluous.

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Once more add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge at 1 min for 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

6. Drying of the Spin Filter

Remove the residual ethanol by final centrifugation for 2 min at 12,000 rpm.

7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 100 µl of the preheated (65°C) Elution Buffer (**Code E**) directly onto the Spin Filter. Incubate on a heating block for 3 min at 65°C (no shaking).

Centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

Note: The DNA can also be eluted with a lower or a higher volume of Elution Buffer (depending on the expected yield of DNA). The minimum volume for the elution is 50 µl.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at +4°C. For a storage time of more than 24 hours it should be kept at -20°C.

Product Information

- Material Safety Data Sheet
 - Validation Report
 - Flow chart
- (Download: <http://www.congen.de/en/company/downloads>)

Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

Distribution and ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com

Safety information

Lysis Buffer

Warning

H319 H412
P280 P273 P305+P351+P338

Proteinase K

Danger

H315 H319 H334 H335
P280 P305+P351+P338

Binding Buffer

Danger

H225 H319 H336
P210 P233 P305+P351+P338

Pre-Wash Buffer

Warning

H302 H412 EUH032
P280 P273 P305+P351+P338

- H225:** Highly flammable liquid and vapour.
H302: Harmful if swallowed.
H315: Causes skin irritation.
H319: Causes serious eye irritation.
H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335: May cause respiratory irritation.
H336: May cause drowsiness or dizziness.
H412: Harmful to aquatic life with long lasting effects.
EUH032: Contact with acids liberates very toxic gas.
P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233: Keep container tightly closed.
P273: Avoid release to the environment.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information we offer a Material Safety Data Sheet. Please contact your distributor or send an e-mail to sales@r-biopharm.de.