# **RIDA<sup>®</sup> QUICK T-2 / HT-2 RQS**

**Art. No.: R5304**

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# RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS

## Brief information

RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS (Art. No. R5304) is an immunochromatographic test for the quantitative determination of T-2 / HT-2 toxin in oats, wheat and corn (maize). All reagents required for the assay are contained in the test kit.

Results are evaluated with the RIDA<sup>®</sup>QUICK SCAN (RQS, Art. No. ZG5005) or the RIDA<sup>®</sup>SMART APP software (Art. No. ZRSAM1000).

Sample preparation: homogenization and extraction

Time requirement: sample preparation (for 10 samples) .....approx. 10 min  
test implementation (incubation time) .....5 min

Limit of detection 50 µg/kg (ppb)

Specificity: The RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS test detects T-2 / HT-2 toxin in oats, wheat and corn (maize) samples. The test has been validated for these matrices.

The specificity of the RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

## Related products

RIDASCREEN<sup>®</sup> T-2 Toxin (R3801)

RIDASCREEN<sup>®</sup>FAST T-2 Toxin (R5302)

RIDASCREEN<sup>®</sup> T-2 / HT-2 Toxin (R3805)

## 1. Intended use

RIDA®QUICK T-2 / HT-2 RQS is a quantitative immunochromatographic test in strip format for the quantitative determination of T-2 / HT-2 toxin in oats, wheat and corn (maize).

## 2. General

T-2 and HT-2 toxin belong to the trichothecene group of mycotoxins and are formed by fungi of the genus *Fusarium*. T-2 and HT-2 toxin are often found in agricultural commodities, although the incidence and the concentrations found show a broad regional variation. Due to their high cytotoxic and immunosuppressive mode of action T-2 / HT-2 toxins are a threat for human and animal health.

## 3. Test principle

The basis of the immunochromatographic assay in test strip format is an antigen-antibody reaction. A specific antibody against T-2 / HT-2 toxin binds the T-2 / HT-2 toxin molecules in the samples.

Results are evaluated with the RIDA®QUICK SCAN (RQS) reader or the RIDA®SMART APP software.

The control band (control line) should be present in all cases in order to prove that the test strip is valid.

## 4. Reagents provided

Each kit contains sufficient materials for 20 determinations.

Each test kit contains:

Component	Cap color	Format	Content
20 x Test strip	-	Ready to use	One determination each (separately packed)
1 x Extraction buffer	transparent	Concentrate 9x	115 ml
1 x Mobile solvent	white	Ready to use	8 ml
1 x Sample buffer	transparent	Ready to use	20 ml
1 x RIDA®SMART APP cover	-	Ready to use	Lot-specific evaluation cover for the use with RIDA®SMART APP

## **5. Materials required but not provided**

### **5.1. Equipment:**

- laboratory mincer / grinder
- balance
- filter paper (Whatman No. 1) or equivalent
- optional: shaker
- optional: centrifuge
- graduated cylinder
- 100 µl, 150 µl, 350 µl and 900 µl micropipettes (e.g. R-Biopharm FP 50 (Art. No. Z0006), FP 100 (Art. No. Z0007) and FP 150 (Art. No. Z0009))
- RIDA<sup>®</sup>QUICK SCAN reader (Art. No. ZG5005) or activated RIDA<sup>®</sup>SMART APP software (Art. No. ZRSAM1000) installed on compatible smartphone

## **6. Warnings and precautions for the users**

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The test strips are sensitive to humidity. Humid test strips may influence the test results negatively; therefore keep the strips away from humidity. This is especially important for already opened test strip packing.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at [www.r-biopharm.com](http://www.r-biopharm.com).

## **7. Storage instructions**

Store the kit at 2 - 8 °C / 36 - 47 °F. DO NOT FREEZE the test strips.

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual components between kits of different lot numbers.

## 8. Sample preparation

For sample preparation a ready to use extraction buffer is needed. For this dilute the extraction buffer (concentrate 9x) 1:9 (1+8) with distilled water. The extraction buffer should be handled under clean working conditions. The ready to use extraction buffer should be stored at 2 - 8 °C (36 - 46 °F) and expires after 8 - 10 weeks.

**Bring test strips and ready to use extraction buffer to room temperature (20 - 25 °C / 68 - 77 °F) before use.**

The samples should be stored in a cool place, protected from light. A representative sample (according to accepted sampling techniques) should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

### 8.1. Sample preparation when performing test evaluating with RIDA®QUICK SCAN reader

#### 8.1.1 Detection limit 50 µg/kg

- weigh 10 g of ground sample into a screw cap tube and add 45 ml of ready to use extraction buffer
- close the tube and shake the sample vigorously for approx. 3 min (manually or with shaker)
- filtrate extract or centrifuge (2 min, 2,000 g)
- mix 150 µl of the clear supernatant with 350 µl mobile solvent
- use 100 µl as sample solution (dilution I) in the test

#### 8.1.2. Detection limit 800 µg/kg

- mix 100 µl of dilution I with 900 µl sample buffer to obtain dilution II
- use 100 µl of dilution II as sample solution in the test

#### **Important:**

**If test strip evaluation is performed with the RIDA®QUICK SCAN reader the obtained result has to be multiplied by the factor 10 to get the actual T-2 / HT-2 toxin concentration of the sample.**

## 8.2. Sample preparation when performing test evaluating with RIDA<sup>®</sup>SMART APP

### 8.2.1 Detection limit 50 µg/kg

- weigh 10 g of ground sample into a screw cap tube and add 45 ml of ready to use extraction buffer
- close the tube and shake the sample vigorously for approx. 3 min (manually or with shaker)
- filtrate extract or centrifuge (2 min, 2000 g)
- mix 150 µl of the clear supernatant with 350 µl mobile solvent
- use 100 µl as sample solution (dilution I) in the test

### 8.2.1. Detection limit 1.000 µg/kg

- mix 100 µl of dilution I with 900 µl sample buffer to obtain dilution II
- use 100 µl of dilution II as sample solution in the test

**If the RIDA<sup>®</sup>SMART APP is used for test strip evaluation please choose the appropriate RIDA<sup>®</sup>SMART APP software application for your sample preparation and matrix (see section 10.1.).**

## 9. Test procedure

- apply 100 µl of sample solution (see section 8.) onto the application area of the test strip
- read the result after 5 min (+ max. 15 sec) incubation with RIDA<sup>®</sup>QUICK SCAN or RIDA<sup>®</sup>SMART APP.

## 10. Evaluation


The right band in the test strip reaction field is a control band (control line) and must appear after each test procedure. If the band is missing, the test result is not valid because of improper test procedure or deterioration of the reagents. Repeat the test with a new strip. If the right band is missing again, please inform your local distributor.

Please note: If the sample is contaminated > 2000 µg/kg visibility of the control band is very weak.

The left band in the reaction field is the test band (test line) and must also appear after each test procedure.

To obtain comparable results, the test strips must always be evaluated after 5 min.

Fig.1: Test strip RIDA<sup>®</sup>QUICK T-2 / HT-2 Toxin RQS

### 10.1. Evaluation with RIDA<sup>®</sup>QUICK SCAN (Art. No. ZG5005)

First, please read the instructions for use of the RIDA<sup>®</sup>QUICK SCAN reader attentively. The test strip evaluation has to be performed according to these guidelines.

Please find a short description regarding the use of RIDA<sup>®</sup>QUICK SCAN in the test kit.

RIDA<sup>®</sup>QUICK SCAN software method: **T2/HT2 Quant**

### 10.2. Evaluation with RIDA<sup>®</sup>SMART APP (Art. No. ZRSAM1000)

First, please read the instructions for use of the RIDA<sup>®</sup>SMART APP software attentively. The test strip evaluation has to be performed according to these guidelines.

Please find a short description of how to use the RIDA<sup>®</sup>SMART APP on our website [www.r-biopharm.de](http://www.r-biopharm.de).

For the test strip evaluation the RIDA<sup>®</sup>SMART APP cover is needed which is included in the test kit. The RIDA<sup>®</sup>SMART APP cover contains lot-specific information and must be used only with the appropriate lot of the R5304 RIDA<sup>®</sup>QUICK T-2 / HT-2 RQS test kit. For further information please check the certificate of analysis enclosed in the kit.

Before starting the test implementation and test evaluation please check if the correct RIDA<sup>®</sup>SMART APP cover is available.

RIDA®SMART APP software applications for the evaluation of RIDA®QUICK T-2 / HT-2 Toxin RQS test strips:

Detection range	Matrices	RIDA®SMART APP application
50 – 1.000 µg/kg	Oats, corn	Oats_Corn 50-1000ppb
1.000 – 10.000 µg/kg	Oats, corn	Oats_Corn 1000-10000ppb
50 – 1.000 µg/kg	Wheat	Wheat 50-1000ppb
1.000 – 10.000 µg/kg	Wheat	Wheat 1000-10000ppb

## 11. Sensitivity

The RIDA®QUICK T-2 / HT-2 RQS test is capable of detecting T-2 / HT-2 toxin contaminations of 50 – 800 µg/kg or 500 – 8,000 µg/kg if test strip evaluation is performed with the RIDA®QUICK SCAN reader. T-2 / HT-2 toxin contaminations of can be detected in the range of  $\geq 50$  - 1000 µg/kg or 1000 - 10000 µg/kg when evaluating with the RIDA®SMART APP.

**For further product information please contact your local distributor or R-Biopharm AG ([sales@r-biopharm.de](mailto:sales@r-biopharm.de)).**

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG