RIDA[®] QUICK Gliadin (ready to swab)

Art. No. R7005

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RIDA[®] QUICK Gliadin (ready to swab)

Brief information

RIDA[®]QUICK Gliadin (ready to swab) (Art. No. R7005) is an immunochromatographic test for the qualitative detection of gliadin / gluten contamination on surfaces. The swab test is used for the hygiene control in production and in laboratories; clean production sites are a prerequisite for gluten-free food.

The identical R5 dip stick RIDA[®]QUICK Gliadin (R7003) has been approved as AOAC-OMA for food samples.

The test kit contains 25 test strips (single packaged) and 25 vials filled with 0.5 ml running buffer for 25 determinations. Results are evaluated visually.

Time requirement: sampling for swab test..... approx. 1 min
 test implementation (incubation time) 5 min

Detection limit: on **surfaces** approx. 1 - 2 µg gliadin / 100 cm²
 (approx. 2 - 4 µg gluten / 100 cm²)

Specificity: The **monoclonal antibody R5** detects the gliadin-fraction from wheat and corresponding prolamins from rye and barley.

Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed / processed food (e.g. maize bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

Related products:

RIDASCREEN[®] Gliadin (Art. No. R7001)

RIDASCREEN[®]FAST Gliadin (Art. No. R7002)

RIDASCREEN[®] Gliadin competitive (Art. No. R7021)

RIDA[®]QUICK Gliadin (Art. No. R7003)

RIDA[®]QUICK Gliadin (single packaged) (Art. No. R7004)

Cocktail (patented) (Art. No. R7006 / R7016)

RIDA[®] Extraction Solution (colorless) (Art. No. R7098)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

SureFood[®] ALLERGEN real time PCR Gluten (Art. No. S3106)

SureFood[®] ALLERGEN QUANT real time PCR Gluten (Art. No. S3206)

1. Intended use

RIDA[®]QUICK Gliadin (ready to swab) (Art. No. R7005) is used as a swab test for the gluten determination on surfaces in the hygiene control. The test has been developed for the detection of low amounts of gluten (contamination). **No** high-dose-hook-effect is observed at high concentrations. However, the red target band may smear at high gluten concentrations.

2. General

Especially airborne cereal dust lead to gliadin contamination in the production. The cleaning of surfaces should be carried out with alcoholic solutions. Thereafter, using swab tests, it is verified if the cleaning of the surfaces was sufficient.

The Codex Alimentarius Commission has stipulated in the „Codex Standard for Foods for Special Dietary Use for Persons Intolerant to Gluten” (CODEX STAN 118-1979) the limit value for gluten-free food at 20 mg/kg gluten. This threshold has also been adopted by many national legislations. The prolamin content (e.g. gliadin) of gluten is generally taken as 50 % (CODEX STAN 118-1979).

The official type I method for gluten determination according to the Codex Alimentarius is an ELISA which uses the R5 antibody (Mendez). This requirement is fulfilled by RIDASCREEN[®] Gliadin test (Art. Nr. R7001). **The RIDA[®]QUICK Gliadin test strips also use the R5 antibody and show a good correlation with the official method, the R5-ELISA RIDASCREEN[®] Gliadin. R-Biopharm AG is the only company, that is allowed to use the R5 antibody for the test strips.**

3. Test principle

The basis of the immunochromatographic test is the monoclonal R5-antibody which is specific for the detection of gliadin from wheat and prolamins from rye and barley. If gliadin is present, a sandwich is formed consisting of immobilized R5 antibody at the target band, gliadin and red latex-labeled R5 antibody. Results are read visually. Generally, the higher the analyte level in the sample, the stronger the red color of the test band will be.

4. Reagents provided

Each kit contains sufficient materials for 25 measurements.

Component	Cap color	Format	Volume
Test strip	-	Ready to use, single packaged	25 piece
Buffer	transparent	Ready to use	25x 0.5 ml
Evaluation card			1 piece

5. Warnings and precautions for the users

Airborne cereal dust may lead to gliadin contamination of the assay. Before carrying out the assay clean surfaces with 40 % ethanol or 2-propanol. Wear gloves to avoid contamination.

The test strips are very sensitive to humidity that could turn the test useless. For this reason keep the strips away from humidity!

This kit may contain further hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

6. Storage instructions

Store the unopened kit at 2 - 8 °C (36 - 46 °F). Do not freeze the kit.

No quality guarantee is accepted after the expiry date on the kit label.

7. Test implementation for swab test

1. Take as many vials as surfaces to be swabbed. Tap the vials onto the table so that the liquid detaches from the lid.
2. Swab the lower end (reaction zone) of a **dry** test strip thoroughly over a sampling area of 10 x 10 cm (wear gloves).
3. Place the test strip vertically into the vial with the arrow end pointing down.
4. Take out the strip after exactly 5 min (+/- 10 s) and read the result using the evaluation card.

8. Results and Sensitivity for Swab Tests

Positive result: two colored bands

The surface is contaminated if two colored bands (the blue control band and the red test band) are visible within the result window. In the case of swabbing test, bands with non-uniform intensity may occur due to an inhomogeneous gluten distribution on the surface or different swabbing procedures.

Concentration > approx. 1 - 2 μg gliadin / 100 cm^2 (approx. 2 - 4 μg gluten/ 100 cm^2)

Negative result: only the blue control band

The surface is gluten-free if no red test band is visible within the result window.

Concentration < approx. 1 - 2 μg gliadin / 100 cm^2 (approx. 2 - 4 μg gluten/ 100 cm^2)

Invalid result: no colored band

If no band is visible within the result window after performing the test, the test is considered invalid.

In general

- The test strip has been developed for the detection of gluten contamination.
- The limit of detection is dependent on the properties of the surface and the kind of contamination.

–A negative result does not necessarily indicate the absence of gluten as the gluten may be inhomogenously distributed or the level of gluten in the product is below the limit of detection.

Recommendation

- For documentation, the upper part of the test strip marked with “Gluten” together with the test bands must be cut off.
- For quality control it is recommended to contaminate the surface with a known amount of gluten and to swab thereafter.
- For screening of raw materials and processed food a sample preparation has to be carried out, for that purpose RIDA[®]QUICK Gliadin (Art. No. R7003) and RIDA[®]QUICK Gliadin (single packaged) (Art. No. R7004) are recommended.
- The ELISA RIDASCREEN[®] Gliadin (Art. No. R7001) should be used for quantification. This test kit is also AOAC-RI and AOAC-OMA (Official Method of Analysis, first action status) validated.

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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- RIDASCREEN® Gliadin (Art. No. R7001)
- RIDASCREEN® FAST Gliadin (Art. No. R7002)
- RIDASCREEN® Gliadin competitive (Art. No. R7021)
- RIDA® QUICK Gliadin (Art. No. R7003)
- RIDA® QUICK Gliadin (single packaged) (Art. No. R7004)
- Cocktail (patented) (Art. No. R7006 / R7016)
- RIDA® Extraction Solution (Art. No. R7098)

Set of 3 processed Gliadin Assay Controls (Art. No. R7012)

SureFood® ALLERGEN real time PCR Gluten (Art. No. S3106)

SureFood® ALLERGEN QUANT real time PCR Gluten (Art. No. S3206)

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are registered trademarks of R-Biopharm AG