User Manual


SureFast[®] PREP Bacteria

Art. No. F1021

100 extractions

Efficient DNA preparation of bacteria

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Table of contents 

Description 7
Principle 7
Kit components (per box) and storage 7
Additionally required equipment and materials 7
Preparations 8
 General 8
 Before each preparation 8
Extraction protocol 8
Product Information 9
Technical Support 9
Safety information.....10

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1x Lysis Buffer (20 ml)	(Code L)	1x Receiver Tubes 2.0 ml, yellow (50 x)	(Code R)
1x Binding Buffer (10 ml)	(Code B)	1x Receiver Tubes 1.5 ml, clear (50 x)	(Code T)
1x Wash Buffer (60 ml)*	(Code W)	1x Spin Filter (50 x)	(Code S)
1x Elution Buffer (10 ml)	(Code E)		

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Description

The kit is intended to be used for the isolation of bacteria - DNA from food (enrichments, avulsions and swabs).

In samples with a very low microbial activity (e.g heat pasteurized egg white powder) a contamination with not viable or not culturable bacteria cells might be observed. In case of other conditions the samples show a fast metabolization of these cells. To assess the process of bacterial growth, it is recommended to compare the samples at the beginning and at the end of the culturing

Principle

1. Preparation of the basic material
2. Lysis at 99°C
3. Setting of optimal binding conditions
4. Binding of the nucleic acids on a Spin Filter
5. Purification of the bound nucleic acids
6. Drying of the Spin Filter
7. Elution of nucleic acids from the Spin Filter

Kit components (per box) and storage

1x Lysis Buffer (20 ml)	(Code L)	1x Receiver Tubes 2.0 ml, yellow (50x)	(Code R)
1x Binding Buffer (10 ml)	(Code B)	1x Receiver Tubes 1.5 ml, clear (50x)	(Code T)
1x Wash Buffer (60 ml)*	(Code W)	1x Spin Filter (50x)	(Code S)
1x Elution Buffer (10 ml)	(Code E)		

* After adding ethanol (purity $\geq 96\%$; not supplied with the kit)

All reagents of the kit should be stored dry and at room temperature (14-25°C).

Additionally required equipment and materials

- suitable equipment for sample comminution and homogenization
- micro balance and spatula for weighing the samples
- reaction tubes free from DNA and DNase 1.5 ml; 2.0 ml
- waterproof pen and tags for labeling the reaction tubes
- unpowdered disposable gloves
- pipettes with filter tips
- Vortex mixer
- Thermomixer/ heating block (up to 99°C)
- micro centrifuge (up to 12,000 rpm)
- ethanol (purity $\geq 96\%$) for preparation of Pre-Wash Buffer and Wash Buffer
- disposal bags or waste bin
- autoclave
- biological safety cabinet

Preparations

General

- Add 42 ml ethanol to the Wash Buffer (**Code W**) and mix thoroughly.

Before each preparation

- Preheating the Elution Buffer (**Code E**) - Transfer the needed amount of Elution Buffer (**Code E**) under calculation of a reserve volume into a reaction tube (not supplied with the kit) and equilibrate to 60°C (The Elution Buffer is necessary for step 7).

Extraction protocol

1. Preparation of the basic material

Transfer 1.0 ml of the enrichment or avulsions under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).

Centrifuge for 5 min at 12,000 rpm.

Discard the fluid and inactivate per sterilization.

For swabs, place them into a 1.5 ml reaction tube (not provided with the kit) and cut the shaft, so that you can close the cap of the extraction tube.

2. Lysis of the basic material

Add 400 µl of Lysis Buffer (**Code L**), mix it briefly and close the reaction tube with caps.

Incubate under continuously shaking for 10 minutes at 99°C.

3. Setting of optimal binding conditions

Centrifuge the sample lysate for 1 min at 12,000 rpm.

Transfer ca. 300 µl of the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit).

Note: Depending on the nature of the sample a smaller amount of liquid may be obtained after centrifugation.

In this case take the actual volume for the following steps or vary the volume of the Lysis Buffer before lysis.

4. Binding of the nucleic acids on a Spin Filter

Add 200 µl Binding Buffer (**Code B**) to the filtrate and mix the sample.

Place a Spin Filter (**Code S**) into a yellow 2.0 Receiver Tube (**Code R**).

Transfer the solution directly onto the filter and incubate at room temperature for 1 min.

Centrifuge for 1 min at 12,000 rpm.

Discard the filtrate and place the Spin Filter back into the yellow Receiver Tube.

5. Purification of the bound nucleic acids

Add 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

Add again 550 µl Wash Buffer (**Code W**) to the Spin Filter and centrifuge for 1 min at 12,000 rpm. Discard the filtrate and place the Spin Filter back into the Receiver Tube.

6. Drying of the Spin Filter

Remove the residual Ethanol by centrifugation for 2 min at 12.000 rpm.

7. Elution of nucleic acids from the Spin Filter

Place the Spin Filter into a clear 1.5 ml Receiver Tube (**Code T**) and add 100 µl of the preheated Elution Buffer (**Code E**) directly onto the Spin Filter.

Incubate for 3 min at room temperature and centrifuge for 1 min at 10,000 rpm. After centrifugation discard the Spin Filter.

The eluted DNA is ready-to-use for the PCR. The DNA can be stored for up to 24 hours at 4°C. For a storage time of more than 24 hours it should be kept at -20°C.

Product Information

- Material Safety Data Sheet
- Flow chart
(Download: <http://www.congen.de/en/company/download>)

Technical Support

For further questions please send an e-mail to info@congen.de.

Safety information

Binding Buffer

Danger

H225-319-336 P210-233-305-351-338

- H225:** Highly flammable liquid and vapour.
H319: Causes serious eye irritation.
H336: May cause drowsiness or dizziness.
P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P233: Keep container tightly closed.
P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

For further information we offer a Material Safety Data Sheet. Please contact your distributor or send an e-mail to info@congen.de.